

Desiccant

Our premium Eurosiv desiccant's are subject to our own comprehensive testing to check: Gas Desorption; Delta-T Results; Initial Moisture Content (should be less than 2.5 %); Moisture Adsorption Capacity (tested to new EN1279 standards); Dust Content and Bulk Density (to check volume).



Eurosiv Standard Bead Desiccant (Box) (1.3 Bead)

Eurosiv Standard	Code	Weight (kg)
Box	D021	25



Eurosiv Extra Mini Bead Desiccant (Box) (0.5 - 0.9mm Bead)

Eurosiv Extra Mini	Code	Weight (kg)
Box	D019	25

NOTE: Not all desiccant boxes of 25kg offer equivalent value for money.

Consider that you fill desiccant by volume and not weight, so look at the size and volume of the box and you'll often find a lower volume than our premium product.

Eurosiv Extra Mini Bead Desiccant (Drum) (0.5 - 0.9mm Bead)



Eurosiv Extra Mini	Code	Weight (kg)
Drum	D024	150

Eurosiv Desiccant Test Kit

Industry Standard test for determining the activity of molecular sieves prior to usage in the process of manufacturing IG units.



Eurosiv Desiccant Test Kit	Code
Delta-T Test Kit	A450

Please order using the product codes and descriptions. If ordering for the first time, please contact: sales@thermosealgroup.com for details of your Area Sales Manager.

Access our Downloads here:

